# AAT088

## Operation Manual

**Breath Alcohol Tester**

Thanks for purchasing our AAT088, a hand-held Digital Alcohol Breath Analyzer. AAT088 gives you a digital reading of your breath alcohol content and your possible concentration of alcohol intoxication.

## How Does Alcohol Testing Work?

When alcohol enters the blood through the alimentary canal (mouth, stomach, and intestines), it is distributed through every part of the body, including the lungs. As the blood passes through the lungs a quantity of alcohol passes through the alveoli (tiny sacs of air in the lungs). This quantity of alcohol that passes through the alveoli is proportionate to the amount of alcohol present in the blood. This period will last for about 15~30min before stable.

By using a breath testing device, it is possible to measure the quantity of alcohol that passes through the alveoli and consequently calculate the person's BAC (Blood Alcohol Concentration). The ratio of breath alcohol to BAC is 2000:1 - in simple terms, every 2000 milliliters of breath will contain the same amount of alcohol as 1 milliliters of blood. (and some other countries is using 2300:1 or 2200:1 or 2100:1)

For this reason, it is necessary for a test subject to blow into the machine for at least 6 seconds so that a sufficient deep air lung sample is obtained. Failure to do so may result in an artificially low reading.

Our device detects human breath and gives a digital reading out of a person's approximate Blood Alcohol Concentration (known from here on as BAC). If user direct blow very high concentration into the sensor, its surface will be surrounded by lots of ethanol particles, which need to wait for 15 min for evaporation.

## ■ Important

To begin, rinse your mouth with fresh water before using the device. Abstain from smoking or drinking for at least 15 minutes. If you perform a test just after drinking, the reading will fluctuate from very high to low reading (i.e. it will drop gradually). It is due to the alcohol residues inside your mouth that effects the test values. This will also happen with aromatic drinks, alcoholic mouthwashes or mouth sprays, medicines or drug, as well as by burping and vomiting.

**Please read the full instructions carefully for proper use of the AAT088 unit. Never drink and drive – the only safe level is zero!**

# ■ Battery Power


AAT088 comes with a built-in 3.7V Li ion battery (300mAH) and please fully charge the tester before use.

**Attention: The reading is invalid if the battery power is low.**

When there is not enough battery power, the “Lo” & “bAt” will flash alternatively. Please charge the tester immediately.

Please connect the provided USB cable and charge the machine. The “LOGO” will turn into red color. The red logo will turn off when the battery is fully charged.

Insert the provided Mouthpiece into the gas inlet located aside of the LCD screen.

# ■ Instruction

## - How to power on and start breathalyzer test

- (1) Press the ⏻ Key to turn on the power and enter into the alcohol mode, “Wait” will turn on.

If the machine has already been turned on, just press ⏻ Key once to activate sensor pre-heating, 【 PH 】 will flash for few times and then start a 19 sec. countdown. When countdown to zero, it is ready to start alcohol test.

- (2) “BLOW” flashes with beep tone and it is ready to take sample. There will has a 5s count down period for auto detect function. You can either do 2.1 Manual Test or 2.2 Active test

### 2.1 Manual Test (if do not drink alcohol or at very low concentration)

Press the ⏻ Key once, “BLOW” lit on. Please start to blow through the mouthpiece steadily and continuously for approximately 5 seconds.

## 2.2 Active Test


Please blow in within the **5s auto detect period**. If alcohol detected, "**BLOW**" lit on. Please blow through the mouthpiece steadily and continuously for approximately 6 seconds.

If there is no alcohol detected within **5s auto detect period**, AAT088 will take about 5s sampling as well and will then show the reading.

While blowing, a long beep will be heard and please keep blowing. When the beep tone stops, sampling finishes and the test result reading will be shown on the display immediately.

### -How to check test record

Under the "**RLC**" mode, press **(i)** key for about 1 seconds to enter into the memory recall mode. Press **(i)** key again to scroll the records. (Max.10 records). Press and hold key for 2 seconds to exit memory recall mode.

In the record mode, the test number and its reading will flash on the LED alternatively. For example, the 9<sup>th</sup> record is **0.58‰ BAC (=0.058% BAC)**, then r09 and 0.58‰ will flash alternatively.

If memories full, the new reading will overwrite the previous one (FIFO). All the records will be erased if user takes off the rechargeable battery.

### -How to power off:

- (1) Auto off, the machine will turn off automatically if no key being press for about 25 seconds.
- (2) Manual, press and hold the **(⏻)** key for 1 second to power off the machine when display show "**RLC**".

### -How to select Unit / Alert setting:

**Note: This is for service center to use, please do not change it unless there are special circumstances.**

- (1) When display show "**RLC**", press "**(i)**" & "**(⏻)**" key at the same

## English

time for 3 seconds and " Unt " will display on the 8-segment display, press the " ⏪ " key to scroll 4 different units.

Unt Unt Unt Unt  
mg/l → ‰ → mg/100mL → ‰

(2) After select the desired unit, either hold the ⓘ for 2 seconds and exit the setting mode, or press the ⓘ once to go to the alert setting mode.

Press the ⏪ key to scroll 4 different alert level, b-2, b-3, b-5, b-8 as following

| Alert value being set |          | Alert level (shown on display) |              |              |              |
|-----------------------|----------|--------------------------------|--------------|--------------|--------------|
|                       |          | b-2                            | b-3          | b-5          | b-8          |
| Unit                  | %BAC     | 0.2‰                           | 0.3‰         | 0.5‰         | 0.8‰         |
|                       | %BAC     | .02% (0.02%)                   | .03% (0.03%) | .05% (0.05%) | .08% (0.08%) |
|                       | mg/l     | 0.1 mg/l                       | 0.15 mg/l    | 0.25 mg/l    | 0.4 mg/l     |
|                       | mg/100mL | 20 mg/100mL                    | 30 mg/100mL  | 50 mg/100mL  | 80 mg/100mL  |

After select the unit and alert level, you can either leave the machine rest for 3 seconds to exit the Unit/Alert setting mode automatically or hold the " ⓘ " key for 2 seconds and exit this setting mode.

### -How to read alert signals

AAT088 alcohol tester is designed to alert user as following:

|              |                                                                                                                                            |
|--------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Low Level    | For measurements between 0~0.08 ‰BAC (0~0.04mg.l), 0 will be displayed and avoid false alarm due to environment or smell in humans' mouth. |
| Safety Level | For measurements between Low Level and Alert Level, the alcohol concentration level shows on the display.                                  |
| Alert Level  | For measurements equal to or above the alert concentration level, the analyzer beeps together "HI" flashes for 10 times                    |

### ■ Caution

- (1) To achieve an accurate reading, user should ensure there is no alcohol residue inside the chamber before starting a new test
- (2) After 4 consecutive tests, please allow the device to rest for 2 min.
- (3) When the product is unused, please store it in a dry environment.
- (4) When low power, the reading is invalid, please recharge the machine immediately.

## English

---

- (5) If the difference between two successive reading is  $>0.1\%$ , user should test again.
- (6) If there is very high concentration tested, user should wait for at least 3 min to preform another test.

## Specification

Display Range: 0/ 0.08~2.00 ‰BAC or 0/ 0.008~0.200 %BAC  
or 0/ 8~200 mg/100mL or 0/ 0.04~1.00 mg/l  
(‰ & mg/l convection factor is about 2.0)

Power Input: 3.7V (rechargeable battery).

Charger: 5V, 2A or above

Display: Digital LED display (3 decimal points)

Record: 10 records

Operating Temperature: 5°C~50°C, for optimum operation

Agency Approval: CE FCC UKCA

Accessories: 10 mouth pieces.

## Precautions

Your blood alcohol concentration (BAC) may change within minutes

Do not use this device to determine whether you can drive or not.

You cannot use the test results in Court.

Your safest bet is not to mix drinking with driving at all. Avoid Binge drinking.

Store it according to specifications. Avoid long term sun exposure and humidity.

## Liability for Proper Function or Damage

Never unscrew the Alcohol Tester. Warranty is not valid if any screws had been loosened or warranty label had been broken; or if the instrument is used in a manner not conforming to its intended use.

Other factors (e.g. environment, smoke & other mouth impurities) may affect the BAC level, please wash your mouth always before testing, wait at least 15 minutes after your last alcohol drink and use the readings for Reference Only.

## Maintenance

Repair/maintenance of the "AAT088 Alcohol Tester" should only be carried out by trained technician from our authorized Service Center. Only authentic "AAT088" spare parts can be used for maintenance.

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